

Video laryngoscope with single-use blade

UED-D



Improve efficiency
Reduce cross-contamination

ZHEJIANG UE MEDICAL CORP.

✉ INFO@UEWORLD.COM

☎ +86 571 28057363

🌐 EN.UEWORLD.COM

📍 8 YOUYI RD, BAITA ECONOMIC
DEVELOPMENT ZONE, XIANJU,
ZHEJIANG 317300, CHINA



UED-M3-S adds on:

Touch screen

- Image mirroring: flips upside-down whenever needed
- Auto shutdown: powers off after 10 mins of inactivity

Wifi ready

- Real-time screen sharing
- Fast access to photos and videos

2 monitor options (UED-M3, UED-M3-S)



UED-D0
newborn



UED-D1
0-2 yrs



UED-D2
2-6 yrs



UED-D3
6-adult



UED-D4
adult, obesity



UED-D5
difficult type



UED-MD1
Miller 1




UED-MD2
Macintosh



UED-MD3
Macintosh

Single-use blades

- Designed to prohibit any reuse



Compatible with UE-M8-S 8" monitor

Grab and go

- Quick set up
- Go-to device for prehospital (EMS) & hospital
- No concerns about maintenance, cleaning, and disinfecting

Whenever you want it

- Lighter
- Longer powered-up
- Faster blades changing (hot-swappable)

Specifications

Monitor		Blades	
Display	3" TFT LCD, 720x480 px	Camera	CMOS
Tilt and rotation	0°-110° tilt, 0°-270° rotate	Angle of view	65°
Memory	32G, take > 400k photos, record ≥16 hrs of video	Light source	high intensity LED
Weight	150g	Weight	40g
Battery life	> 3.5 hrs on a full charge		
Charging time	< 3 hrs		